

# WHICH TABLE IS RIGHT FOR YOU?

## Urology C-Arm Table - 800

Configure to your exact specifications



### SPECIFICATIONS:

Dimensions: 84.5" l x 25.5" w (214.6 x 64.8 cm) with OR accessory rails  
Tabletop: 53" l x 21" w (134.6 x 53.3 cm); 84" l (213.4 cm) with extension  
Lightweight Extension (for patient transfer only, not for imaging): 26.5" l x 21" w (67.3 x 53.3 cm)  
Radiolucent Extension: 31" l x 25.5" w (78.7 x 64.8 cm) with OR accessory rails  
Radiolucent Area: 54" l x 17.5" w (137.2 x 44.5 cm)  
Tabletop Material: Three Fibersin™ sections – each with individual 2" mattress  
Table Sections: Four Sections: Adjustable Fowler back 0° to 80°; abdominal section with cutout; radiolucent removable extension and a lightweight non-imaging transfer extension for transferring patients on/off the table  
Mattress: Individual 2" thick (5 cm), one for each section

### Motions:

Height Adjustable: 34" to 43.5" (86.4 to 110.5 cm)  
Tabletop X Motion (head-to-toe): 20" (51 cm)  
Tabletop Y Motion (side-to-side): 10" (25.4 cm)  
Trendelenburg: 0° to 20°  
Reverse Trendelenburg: 0° to 20°  
Lateral Roll: 0° to 20°  
Fowler Back: Adjusts from 0° to 80° with touch handle lock/unlock

### Controls:

Hand Control: activates all motions  
Foot Control: activates all motions  
Main Control Panel: LCD screen displays actual tabletop position (all axes) in real time, and provides menus and prompts for other functions  
Joystick Controller: for motorized X-Y positioning

X-Ray Film Tray: Removable film cassette tray for 14" x 17" film with cassette mask for 10" x 12" film. Film cassette tray can be removed from either side of the table.

X-Ray Grid: Removable X-Ray grid, 103 line, 8:1 ratio

### Wheels:

Foot End – 3" (7.62 cm) locking swivel casters  
Head End – 5" (12.7 cm) swivel casters with central locking  
Accessory Rails: Standard OR accessory rails 1.125" x .375" (2.86 x .95 cm) on table (rails also on table imaging extension)  
Patient Restraints: Two body straps, two ankle straps  
Finish: Powder coat  
Patient Capacity: 400 lb (181.4 kg); weight tested to four times the patient load rating.  
Power: 115 VAC or 230 VAC and Battery. Battery automatically charges when table is plugged in. An additional battery can be charged with optional wall-mounted charger.  
Warranty: One year  
Certifications: UL-2601-1®, CAN/CSA® C22.2 No.: 601-1-M90, EN 60601-1, IEC 60601-1-2:2001-09

**058-800** Table, Urology C-Arm - 800, 115 VAC  
**058-805** Table, Urology C-Arm - 800, 230 VAC

### Optional:

**056-851** Drain Bag Frame  
**058-860** Drain Pain, Stainless Steel  
**056-853** Crutch, Knee, Simplicity, (pair)  
**058-862** Stirrups, PAL with boot pad (pair)  
**058-865** Arm Board, Standard (pair)  
**058-863** Arm Board, Dura Board with pad (pair)  
**056-850** IV Pole, rail mounted  
**056-862** Battery and Wall Charging Unit, 115 VAC  
**056-863** Battery and Wall Charging Unit, 230 VAC  
**056-463** Grid, X-ray, 103-line 8:1 ratio, 18" x 18"  
**058-866** Cassette Tray Holder, 11" x 14"  
**058-867** Extension Arm, Controller

## Brachytherapy C-Arm Table - 810

For image-guided seed implantation



### SPECIFICATIONS:

Dimensions: 84.5" l x 25.5" w (214.6 x 64.8 cm) with OR accessory rails  
Tabletop: 53" l x 21" w (134.6 x 53.3 cm); 84" l (213.4 cm) with extension  
Lightweight Extension (For patient transfer only, not for imaging.): 26.5" l x 21" w (67.3 x 53.3 cm)  
Radiolucent Extension: 31" l x 25.5" w (78.7 x 64.8 cm) with OR accessory rails  
Radiolucent Area: 54" l x 17.5" w (137.2 x 44.5 cm)  
Tabletop Material: Three Fibersin™ sections – each with individual 2" mattress  
Table Sections: Four Sections: Adjustable Fowler back 0° to 80°; abdominal section with cutout; radiolucent removable extension and a lightweight non-imaging transfer extension for transferring patients on/off the table  
Mattress: Individual 2" thick (5 cm), one for each section

### Motions:

Height Adjustable: 34" to 43.5" (86.4 to 110.5 cm)  
Tabletop X Motion (head-to-toe): 20" (51 cm)  
Tabletop Y Motion (side-to-side): 10" (25.4 cm)  
Trendelenburg: 0° to 20°  
Reverse Trendelenburg: 0° to 20°  
Lateral Roll: 0° to 20°  
Fowler Back: Adjusts from 0° to 80° with touch handle lock/unlock

### Controls:

Hand Control: activates all motions  
Foot Control: activates all motions  
Main Control Panel: LCD screen displays actual tabletop position (all axes) in real time, and provides menus and prompts for other functions  
Joystick Controller: for motorized X-Y positioning

X-Ray Film Tray: Removable film cassette tray for 14" x 17" film with cassette mask for 10" x 12" film. Film cassette tray can be removed from either side of the table.

X-Ray Grid: Removable X-Ray grid, 103 line, 8:1 ratio

### Wheels:

Foot End – 3" (7.62 cm) locking swivel casters  
Head End – 5" (12.7 cm) swivel casters with central locking  
Accessory Rails: Standard OR accessory rails 1.125" x .375" (2.86 x .95 cm) on table (rails also on table imaging extension)  
Patient Restraints: Two body straps, two ankle straps  
Finish: Powder coat  
Patient Capacity: 400 lb (181.4 kg); weight tested to four times the patient load rating.  
Power: 115 VAC or 230 VAC and Battery. Battery automatically charges when table is plugged in. An additional battery can be charged with optional wall-mounted charger.  
Warranty: One year  
Certifications: UL-2601-1®, CAN/CSA® C22.2 No.: 601-1-M90, EN 60601-1, IEC 60601-1-2:2001-09

**058-810** Table, Brachytherapy - 810, 115 VAC  
**058-815** Table, Brachytherapy - 810, 230 VAC

### Optional:

**056-853** Crutch, Knee, Simplicity, (pair)  
**058-862** Stirrups, PAL with boot pad (pair)  
**058-865** Arm Board, Standard (pair)  
**058-863** Arm Board, Dura Board with pad (pair)  
**056-850** IV Pole, rail mounted  
**056-862** Battery and Wall Charging Unit, 115 VAC  
**056-863** Battery and Wall Charging Unit, 230 VAC  
**056-463** Grid, X-ray, 103-line 8:1 ratio, 18" x 18"  
**058-866** Cassette Tray Holder, 11" x 14"  
**058-867** Extension Arm, Controller

**NEW**

## 3D Brachytherapy C-Arm Table - 820

Designed to function with a 3D C-Arm



### SPECIFICATIONS:

Dimensions: 103.75" l x 26" w (264 x 66 cm) with OR accessory rails  
Tabletop: 103.75" l x 24" w (264 x 61 cm) includes head and leg extension  
Radiolucent Area:  
3D Imaging Area: 22.5" l x 24" w (57.1 x 61 cm)  
Radiolucent Extension (Fibersin™, non-3D): 29" l x 22.5" w (73.6 x 57.1 cm)  
Total Radiolucent Area (Combined): 51.5" l (131 cm)  
Head-End Extension (non viewable): 18.2" l x 24" w (46 x 51 cm)  
Base Clearance to Floor: 1.5" h (4 cm)  
Tabletop Material: Carbon fiber  
Mattress: Vinyl covered 2" thick (5 cm)

### Motions:

Height Adjustable: 34" h to 44" h (86 to 112 cm)  
Tabletop X Motion (head to toe): 10" (25 cm)  
Tabletop Y Motion (side to side): 10" (25 cm)  
Trendelenburg: 0° to 20°  
Reverse Trendelenburg: 0° to 20°  
Lateral Roll: 0° to 20°

### Controls:

Hand Control: activates all motions  
Foot Control: activates all motions  
Main Control Panel: LCD screen displays actual tabletop position (all axes) in real time, and provides menus and prompts for other functions.  
Joystick Control: for motorized X-Y positioning

### Wheels:

Head End: 3" (7.6 cm) swivel casters, integral locking system on base  
Foot End: 5" (12.7 cm) swivel casters with central locking  
Accessory Rails: Standard OR accessory rails, 1.125" x .375" (2.86 x .95 cm). Head-End 30" l (75 cm); Foot-End 7.75" l (20 cm)  
Patient Capacity: 500 lb (227 kg); weight tested to four times the patient load rating.  
Power: 115 VAC or 230 VAC and Battery. Battery automatically charges when table is plugged in. An additional battery can be charged with optional wall-mounted charger.  
Weight: 542 lb (246 kg)  
Shipping Weight: 684 lb (310 kg)  
Warranty: One year  
Certifications: ETL and cETL listed to UL-60601-1, CAN/CSA® C22.2 No.: 601-1-M90, EN 60601-1 and IEC- 60601-1-2:2001-09

**058-820** Table, 3D Brachytherapy C-Arm, 115 VAC  
**058-825** Table, 3D Brachytherapy C-Arm, 230 VAC

### Optional:

**056-862** Battery and Wall Charging Unit, 115 V  
**056-863** Battery and Wall Charging Unit, 230V  
**056-851** Drain Bag Frame  
**056-850** IV Pole, rail mounted  
**056-853** Crutch, Knee, Simplicity (pair)  
**058-862** Stirrups, PAL with boot pad (pair)  
**058-863** Arm Board, Dura Board with pad (pair)  
**058-865** Arm Board, Standard (pair)  
**058-867** Extension Arm, Controller